

The Advantage of Combining Resin Containing Medium with Saponin Supplemented Lytic Medium



Blood cultures represent one of the most important specimens submitted to the Microbiology Laboratory. The mortality due to sepsis remains high despite appropriate antimicrobial therapy. Therefore, it is important, that laboratories maintain rapid and reliable methods for detection of microorganisms in blood.¹ Detection and identification of anaerobic bacteria in blood cultures is a well-recognized challenge in clinical microbiology. This may be explained by the fact that these microorganisms are typically fastidious, slow growing and difficult to culture.² Furthermore, the anaerobic bottle may continue to provide clinically relevant information because other studies have demonstrated that 4.2 to 9.4% of significant isolates from blood were anaerobes.¹ However, the developments of BD BACTEC™ Automated Blood Culture System and anaerobic blood culture bottles have improved the detection of these microorganisms.²

BACTEC™ Lytic/10 Anaerobic/F culture bottles (pre-reduced enriched Soybean-Casein Digest broth with CO₂) are for anaerobic blood cultures. The principle used is with the **BACTEC** fluorescent series instruments for the qualitative culture and recovery of anaerobic microorganisms from blood. The Lytic anaerobic/F bottle contains saponin as a lysing agent to lyse blood cells and release intracellular microorganisms to facilitate early pathogen detection. Hollick G.E. et al. demonstrated that the addition of saponin to anaerobic media has been shown to enhance the recovery of both anaerobic and facultative anaerobic bacteria.³




Figure 1. Structure of the Saponin

Lysis of blood cells which allows release of viable intracellular organisms or growth factors contained in the cells

In an early study, Rohner and colleagues have proved that combination of Plus aerobic/F and Lytic anaerobic/F was a valuable pair of blood culture media. Plus aerobic/F performs better for patients under antibiotic treatment, due to the antimicrobial-neutralizing effect of resins. For patients without antibiotic therapy, more organisms could be isolated from Lytic anaerobic/F due to cell lysis. They also reported more *S. aureus* isolates ($P = 0.05$), *Pseudomonas spp.* ($P < 0.0001$), other gram-negative bacteria ($P = 0.004$) and yeasts ($P < 0.0001$) were isolated from Plus aerobic/F medium, but more streptococci ($P < 0.0001$), *E. coli* ($P = 0.02$) strains and anaerobes ($P < 0.0001$) were detected with Lytic anaerobic/F medium. Lytic anaerobic/F bottles were significantly ($P < 0.0001$) more often positive at least 6 hours before Plus aerobic/F bottles ($n = 112$ vs. 52, respectively).⁴

A recent study (Rocchetti A. et al, 2016) have evaluated the positivity rate and time to detection of three BACTEC media in an emergency department of Alessandria General Hospital in Italy: Plus Aerobic/F, Plus Anaerobic/F and Anaerobic Lytic/F. During an 8- month period between July 2014 and March 2015, 1503 blood culture sets were collected from 688 adult patients (> 18 years of age).⁵

Table 1
Positivity rate and time to detection.

Culture media	N. ^a	N.Tot ^b	% ^c	Average TTP ^d	Median	Sd ^e
Single media						
Anaerobic media type	177	255	69.4%	21.5	15.4	17.2
Aerobic media type	195	255	76.5%	17.6	14.5	12.3
Lytic media type	196	255	76.9%	16.3	12.0	14.4
Paired media (at least one of the two media type is positive)						
Lytic + aerobic media type	241	255	94.5%	15.7	12.3	11.3
Anaerobic + aerobic media type	224	255	87.8%	18.4	14.1	14.8
Paired media (when the two media type are positive together)						
Lytic + aerobic media type	150	255	58.8%	12.3	11.2	5.9
Anaerobic + aerobic media type	148	255	58.0%	13.6	13.0	4.3

^a Number of positive vials (for single media) and positive sets (for paired media).

^b Total Number of positive sets; on a total of 1503 sets for blood cultures 255 resulted positive for microbiological growth. Sets were composed of anaerobic medium, Aerobic medium and lytic medium. A set was considered positive if at least one of the vials was positive.

^c Percent of positive media on a total of 255 positive sets for blood culture.

^d Average TTP (ttp): time from incubation to time to positivity.

^e Standard deviation.

- Positivity rate (PR) was similar for Anaerobic Lytic/F bottle (76.9%) and Plus Aerobic/F bottle (76.5%) but better if compared with Plus Anaerobic/F bottle (69.4%).
- Median time to detection (TTD) was optimum with the Anaerobic Lytic/F bottle (12 h) versus Plus Aerobic/F bottle (14.5 h) and Plus Anaerobic/F bottle (15.4 h).
- The PR and median TTD for the combination of Plus Aerobic/F with Anaerobic Lytic/F (94.5% and 12.3h, respectively) was significantly better than Plus Aerobic/F with Plus Anaerobic/F (87.8 and 14.1h).



Table 2

Recovery rate and percentage of microorganisms identified in each media. Percentages in brackets indicate percentage of total number for that isolate (N. Tot).

Organisms class	Organisms	N.Tot	Lytic_NP	Anaerobic_NP	Aerobic_NP	
Anaerobic bacteria	<i>Bacteroides fragilis</i>	3	2 (66.7%)	1 (33.3%)		
	<i>Bacteroides stercoris</i>	2	2 (100%)			
	<i>Bacteroides caccae</i>	1	1 (100%)			
Gram-negative bacilli	<i>Escherichia coli</i>	96	77 (80.2%)	62 (64.6%)	69 (71.9%)	
	<i>Klebsiella pneumoniae ssp pneumoniae</i>	14	12 (85.7%)	10 (71.4%)	13 (92.9%)	
	<i>Proteus mirabilis</i>	13	12 (92.3%)	10 (76.9%)	8 (61.5%)	
	<i>Providencia stuartii</i>	9	9 (100%)	9 (100%)	8 (88.9%)	
	<i>Pseudomonas aeruginosa</i>	7	1 (14.3%)	2 (28.6%)	7 (100%)	
	<i>Acinetobacter baumannii</i>	3	1 (33.3%)	1 (33.3%)	3 (100%)	
	<i>Klebsiella oxytoca</i>	3	3 (100%)	3 (100%)	3 (100%)	
	<i>Enterobacter aerogenes</i>	2	2 (100%)	2 (100%)	2 (100%)	
	<i>Enterobacter cloacae</i>	2	1 (50%)	1 (50%)	1 (50%)	
	<i>Salmonella enterica</i>	2	2 (100%)	2 (100%)	2 (100%)	
	<i>Aeromonas salmonicida</i>	1		1 (100%)		
	<i>Aeromonas hydrophila</i>	1		1 (100%)	1 (100%)	
	<i>Brevundimonas vesicularis</i>	1			1 (100%)	
	Gram-positive cocci	<i>Staphylococcus aureus</i>	41	36 (87.8%)	37 (90.2%)	36 (87.8%)
		<i>Staphylococcus epidermidis</i>	15	11 (73.3%)	10 (66.7%)	12 (80%)
		<i>Enterococcus faecalis</i>	13	12 (92.3%)	12 (92.3%)	12 (92.3%)
<i>Streptococcus pneumoniae</i>		8	6 (75%)	8 (100%)	8 (100%)	
<i>Staphylococcus hominis ssp. hominis</i>		6	2 (33.3%)	2 (33.3%)	4 (66.7%)	
<i>Streptococcus gordonii</i>		4	2 (50%)	3 (75%)		
<i>Streptococcus mitis</i>		3	1 (33.3%)	2 (66.7%)		
<i>Allotococcus otitidis</i>		2	2 (100%)			
<i>Gemella morbillorum</i>		2		1 (50%)	2 (100%)	
<i>Peptoniphilus asaccharolyticus</i>		2	2 (100%)	1 (50%)		
<i>Streptococcus anginosus</i>		2	2 (100%)	2 (100%)	2 (100%)	
<i>Staphylococcus capitis</i>		1	1 (100%)	1 (100%)	1 (100%)	
<i>Streptococcus sanguinis</i>		1			1 (100%)	
Other Gram positive bacteria/yeasts		<i>Candida parapsilosis</i>	5	2 (40%)	2 (40%)	5 (100%)
	<i>Listeria monocytogenes</i>	2	2 (100%)	2 (100%)	2 (100%)	
	<i>Candida albicans</i>	2	1 (50%)		2 (100%)	
	<i>Candida glabrata</i>	1	1 (100%)		1 (100%)	
	Recovery rate (for each media) ^a		76.3%	70.0%	77.0%	

^a On a total of 270 isolates.

- Improved recovery of strict anaerobes, facultative anaerobes (gram-negative rods) and yeasts was observed with the Anaerobic Lytic/ F bottle vs. the Plus Anaerobic/ F bottle
- Recovery of gram-positive cocci was observed to be equivalent in the Anaerobic Lytic/ F bottle and the Plus Anaerobic bottle

Table 3

Positivity rate and time to detection for patients receiving antibiotic therapy.

Culture media	N. ^a	N.Tot ^b	% ^c	Average		
				TTP ^d	Median	Sd ^e
Single media						
Anaerobic media type	28	40	70%	27.6	23.3	18.3
Aerobic media type	29	40	72.5%	24.1	18.4	18.5
Lytic media type	30	40	75%	18.0	15.1	8.1
Paired media (at least one of the two media type is positive)						
lytic + aerobic media type	37	40	92.5%	19.7	17.2	14.7
Anaerobic + aerobic media type	36	40	90%	24.6	18.5	19.4
Paired media (when the two media type are positive together)						
No data available						

^a Number of positive vials (for single media) and positive sets (for paired media).^b Total number of positive sets; on a total of 1503 sets for blood cultures 255 resulted positive for microbiological growth. Sets were composed of anaerobic medium, aerobic medium and lytic medium. A set was considered positive if at least one of the vials was positive.^c Percent of positive media on a total of 255 positive sets for blood culture.^d Average TTP (ttp): time from incubation to time to positivity.^e Standard deviation.**Conclusion:**

- The use of resin containing medium with a saponin-supplemented lytic media further improve overall recovery and earlier detection of bacteria and yeasts in the blood of septic patients.
- Differences in performance were evident when comparing Anaerobic Lytic/F media vs. Plus Anaerobic/F media: significantly better recovery of strict anaerobes, strict aerobic & facultative anaerobic gram-negative rods and yeasts was observed with the Anaerobic Lytic/F media.
- Median TTD was significantly decreased with Anaerobic Lytic/F (2.6 h faster for all patients and 3.3 h faster for patients receiving antibiotics). This is because the lysing ability of the Anaerobic Lytic/F media, which reduces background signal due to host cell metabolic activity.

References

1. Hollick G.E. et al. Clinical Comparison of the BACTEC 9000 Standard Anaerobic/F and Lytic/ F Blood Culture Media. *Diagn Microbiol. Infect Dis* 1996;24:191-196.
2. Almuhayawi M. et al. The Performance of the Four Anaerobic Blood Culture Bottles BacT/ALERT-FN, -FN Plus, BACTE-Plus and –Lytic in Detection of Anaerobic Bacteria and Identification by Direct MALDI-TOF MS. *PLoS ONE* 10(11): e0142398. doi:10.1371/journal.pone.0142398.
3. Package Insert: BACTEC™ Lytic/10 Anaerobic/F Culture Vials - Soybean-Casein Digest Broth in a Plastic Vial.
4. Rohner P. et al. Advantage of Combining of Resin with Lytic BACTEC Blood Culture Media. *J Clin. Microbiol.* 1997;35 (10):2634-2638.
5. Rochetti A. et al. Prospective Study of the Clinical Performance of Three BACTEC Media in a Modern Emergency Department: Plus Aerobic/F, Plus Anaerobic/F, and Anaerobic Lytic/F. *Journal of Microbiological Method* 130 (2016) 129-132.



BMS DIAGNOSTICS (M) SDN BHD
(Company Reg. No.: 485573V, GST Reg. No.:001836843008)
19, Jalan 4/62A, Bandar Menjalara, Kepong, 52200 Kuala Lumpur, Malaysia.
Website: www.bmsd.com.my Email: info@bmsd.com.my
Tel: +603- 6272 0236 Fax: +603- 6277 0750

